

ACME™ HEATING PAD



INTRODUCTION

Thank you for purchasing an ACME Heating Pad. With proper care and safe handling, this product will provide relief to you and your loved ones for quite some time.

This pad is designed with the consumer in mind. It has state-of-the-art efficiency and safety features, making it safe and easy-to-use.

PLEASE READ THESE INSTRUCTIONS CAREFULLY. Pay special attention to all warnings and precautionary statements. Failure to adhere to these instructions could result in damage to the product, voiding the warranty; failure to understand the warnings and precautions could make the device unsafe to use, resulting in property damage, serious injury, or death.

GETTING STARTED

In the box, you'll find:

- Instruction booklet
- Warranty information
- Rubber heating mat
- Inner cloth cover
- Outer cloth cover

Caution: Do not pull on the electrical cord to remove the heating pad from the box.

SAFETY WARNINGS AND PRECAUTIONS:

Consult a physician before using this device for prolonged periods of time.

DO NOT use this device on infants, sleeping persons, comatose patients, or people with poor circulation.

DO NOT use in oxygen-rich environments.

NEVER lay or sit on top of the heating pad.

NEVER place the rubber heating mat directly on skin. ALWAYS have both inner and outer cloth covers attached while the heating pad is in use.

ALWAYS before turning on the heating pad, check the rubber mat and electrical cord for tears or holes. DO NOT use this device if the rubber mat or the electrical cord have been folded, crimped, or otherwise compromised. DO NOT use the device if it is wet. **Failure to acknowledge damage to the device could result in blisters, burns, fire, or electrical shock.**

OPERATING THE DEVICE

When you have identified an injured area that needs heating, perform the following:

1. Attach the inner and outer covers to the rubber mat.
2. Plug the electric cord into a 110-volt outlet.
3. Use the tie strings to hold the covered pad on top of the area to be treated. **NEVER use pins or other metal fasteners to secure the heating pad in place.**
4. The device has three (3) settings: Low, Medium, and High.
5. Use the pad for 10-minutes on each setting, starting on "Low" and ending on "High."
6. After 30-minutes of total heating, remove the pad from the area of treatment.
7. Allow the pad to cool down completely before treating another area or putting the device back in the box.

For help with troubleshooting issues, or if you are seeking further clarification of this product's instructions, safety precautions, and warranty information, call 1-800-HOT-PADS during regular business hours to speak with a representative who can assist you.

Visit our website:

www.ACME.com

For instructional videos and customer forums, go to:

www.ACME.com/products/heatingpad/

"At ACME, safety goes in before the heat goes on!"